


# Patient Instructions

## STEP 1 Sample Deposit

- A** Lift toilet seat and position sample collection paper across the rim of the toilet bowl. Secure adhesive tabs to the sides of toilet rim. Lower the toilet seat.




- B** Make bowel movement onto collection paper.




**Attention!**

Fill out personal information on attached label on tubes.




## STEP 2 Sample Collection


- A** Unscrew the purple cap from the sample collection tube. **DO NOT POUR OUT THE LIQUID.**



- B** Poke spiral applicator into stool at 6 different sites. Use only enough fecal material to cover the tip of the applicator. **DO NOT CLUMP, SCOOP, OR FILL THE TUBE.**




- C** Screw the applicator back into the tube and secure tightly.




## STEP 3 Sample Return

- A** Complete the address return envelope to your doctor or laboratory.



- B** Insert sample collection tube into specimen pouch and seal. Insert specimen pouch into return envelope and seal.



- C** Return the sample packet immediately to the clinic or laboratory by mail or in person.

